



OSTEONECTINA



**noraylims**

Integration brand

## What's ServiLIMS?

Software for the management of scientific and technical services of universities and research centres.

ServiLIMS enables both management and web-communication among researchers, service providers and centre administration. ServiLIMS can be adapted to the requirement of the customer by means of the creation and addition of new features.

## WEB PART

### 1 Public Part

- Home page
- Offered services
- Presentation and contact



### 2 Private - Login required

- Service request forms
- Request traceability
- See the results

## LOCAL PART

### 3 Specific modules

- Request Management
- Results Management
- Invoicing module
- Material storage

### General modules

- Statistics and reports
- User management
- Quality-Change control
- Configuration



## ServiLIMS has 3 main parts:

1. Public Web Part: a web portal to show the group description, offered services, price list, etc.
2. Private Web Part: logged users will be able to access here to fill request forms, to trace request state and to see result files.
3. Local Part: software for service managers, to manage requests, register obtained results and create invoice orders.

## Advantages

---

- NorayBio is the first Spanish bioinformatics company and we are specialized in the design and development of software for Biosciences, with high knowledge in LIMS.
- NorayLIMS enables to create a complete central management system as a base to support, in the future, the complete integration of new processes and activities or software.
- NorayBio offers the possibility of developing tailored projects so that the customer can have his own completely customized software.
- Sample traceability maintenance.
- Central and standard information management.
- User-friendly software.

## Technical features

---

NorayLIMS works locally installed or as client/server application.

### Server requirements:

#### • **Minimum:**

- Windows: NT, 2000 or 2003
- Intel Pentium III or AMD Athlon 800MHz
- 512 MB RAM
- 500 MB of free space in the Hard Disk
- Net Ethernet 100 MB

#### • **Recommended:**

- Windows: 200 or 2003 Server
- Intel Pentium IV or AMD XP 1.8GHz
- 1 GB RAM
- 2 GB of free space in the Hard Disk
- Net Ethernet 1 GB

### Workstation requirements:

#### • **Minimum:**

- Windows: NT, 2000 or XP
- Intel Pentium 500MHz
- 128 MB RAM
- 100 MB of free space in the Hard Disk

#### • **Recommended:**

- Windows: 2000 or XP
- Intel Pentium III 1GHz
- 512 MB RAM
- 150 MB of free space in the Hard Disk
- Net Ethernet 1 GB
- Screen resolution: 1024 x 768



### Contact information

NorayBio  
Parque Tecnológico de Bizkaia, 801A  
48160 Derio (Bizkaia) SPAIN  
T.:+34 944036998  
F.:+34 944036999  
info@noraybio.com  
www.noraybio.com